DE DEUTSCH

1. Produktbeschreibung

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EN ENGLISH

1. Product description

Slident zirconium dioxide milling blanks are made of yttrium-stabilized zirconium dioxide (γ -TZP) for dental application of type II, class 4 (Zirkon BioStar HT Smile) or class 5 (Zirkon BioStar Ultra / Zirkon BioStar) in accordance with DIN EN ISO 6872, and they fulfill their material-specific requirements after specified final sintering.

Depending on material grade, Slident zirconium dioxide exhibits different levels of translucency. To enable the dentist to provide each patient with customized and aesthetic restorations, Slident zirconium dioxide milling blanks are available unglazed and precolored.

2. Intended purpose

Slident zirconium dioxide milling blanks are intended for the fabrication of fixed restorations for long-term use.

3. Indication

Zirkon BioStar (colour)

For fabrication of reduced crowns, reduced bridges***, hybrid abutments, bars and cantilever bridges*** as anterior and posterior restorations.

Zirkon BioStar Ultra (colour, ML)

For fabrication of veneers, inlays, onlays, reduced crowns, monolithic crowns, hybrid abutment crowns, monolithic and reduced bridges*** and cantilever bridges*** as anterior and posterior restorations.

Zirkon BioStar HT Smile (colour, ML)

For fabrication of veneers, inlays, onlays, reduced crowns, monolithic crowns, hybrid abutment crowns, monolithic and reduced bridges (up to 3-unit) as anterior and posterior restorations.

*** Up to two adjacent pontics.

*** In Canada the indication for bridges is limited to a maximum of six units with up to two pontics.

*** Cantilever pontics must not be dimensioned longer than 2/3 of the load-bearing abutment crown from mesial to distal.

4. Contraindications

Parafunction, insufficient space, unsuitable preparation, insufficient dental hard tissue, intolerance to components and inadequate oral hygiene.

5. Intended users

Slident zirconium dioxide milling blanks must only be used by dentists and dental technicians in compliance with the procedures in the Instructions for Use.

6. Intended patient group

Fixed restorations using Slident zirconium dioxide milling blanks are suitable for the permanent dentition in adult patients of any gender and nationality.

7. Handling and storage

Before using the material for the first time, check that the packaging and the blank itself are intact. Check whether the contents of the packaging correspond to the declaration on the label. Never use damaged material. Store the blanks only in the original packaging and in a cool, dry place. Avoid vibration, contamination and all contact with fluids.

8. Instructions for use in the laboratory

8.1 Processing / construction
Slident zirconium dioxide is a sensitive, high-performance ceramic and should be processed with special care, including in the partially sintered state!

The following construction parameters must always be considered when fabricating zirconium dioxide constructions:

Indication	Zirkon BioStar (colour), Zirkon BioStar Ultra (colour, ML)		Zirkon BioStar HT Smile (colour, ML)	
	Minimum wall thickness [mm]	Connector cross section [mm ²]	Minimum wall thickness [mm]	Connector cross section [mm ²]
Single crown	okkusal	0,5	0,5	–
	okkusal	0,5	0,5	–
Telescope	okkusal	0,7	0,7	–
	okkusal	0,5	0,5	–
Anterior Bridge 3 units	okkusal	0,5	0,8	> 10
	okkusal	0,5	0,8	> 12
Posterior Bridge 3 units	okkusal	0,7	0,8	–
	okkusal	0,7	0,8	> 12
Anterior Bridge 4 or more units	okkusal	0,6	0,8	> 9
	okkusal	0,8	0,8	> 12
Posterior Bridge 4 or more units	okkusal	0,7	0,8	–
	okkusal	0,8	0,8	> 12
Cantilever Bridge with 1 Pontic	okkusal	0,7	0,8	> 12
	okkusal	0,7	0,8	> 12

Please note that the connector cross-section may have to be of larger dimensions, depending on the construction. For example, in the case of long-span posterior bridges (made of Zirkon BioStar or Zirkon BioStar Ultra), the connector cross-section between two pontics should be increased to at least 20 mm² if possible. The aim is to achieve an oval connector cross-section; the height of the connector is decisive for stability. The design should avoid sharp edges and acute angles. Frameworks for ceramic veneering should be designed in such a way that they support the veneering ceramic in the area of the cusps and embrasures.

Nesting recommendations:
For the design and positioning of the construction in the multilayer blank (Zirkon BioStar HT Smile ML and Zirkon BioStar Ultra ML), the individual layer thicknesses can be taken from the following nesting table:

Blank height [mm]	Zirkon BioStar Ultra ML				Zirkon BioStar HT Smile ML			
	Layer 1+2: Incisal (mm / %)	Layer 3: Intermediate layer (mm / %)	Layer 4: Circumferential layer (mm / %)	Layer 5: Dentin (mm / %)	Layer 1+2: Incisal (mm / %)	Layer 3: Intermediate layer (mm / %)	Layer 4: Circumferential layer (mm / %)	Layer 5: Dentin (mm / %)
14 mm	3,5 / 24,9	2,1 / 15	2,1 / 15	6,3 / 45,1	3,5 / 24,9	2,1 / 15	2,1 / 15	6,3 / 45,1
18 mm	3,5 / 19,4	2,1 / 11,7	2,1 / 11,7	10,3 / 57,2	3,5 / 19,4	2,1 / 11,7	2,1 / 11,7	10,3 / 57,2
22 mm	3,5 / 15,9	2,1 / 9,6	2,1 / 9,6	14,3 / 64,9	3,5 / 15,9	2,1 / 9,6	2,1 / 9,6	14,3 / 64,9
25 mm	3,5 / 14	2,1 / 8,4	3,3 / 13,2	16,1 / 64,4	3,5 / 14	2,1 / 8,4	3,3 / 13,2	16,1 / 64,4

* measured according to DIN EN ISO 6872

11. Possible side effects and interactions
No known side effects or interactions.

12. Disposal
In compliance with local regulations. Non-contaminated and completely emptied packaging can be recycled.

Please note: Observe the information in the latest version of the safety data sheet.

13. Reporting incidents
Any serious incident that occurs in connection with the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

Please note: The safety and clinical performance summary report can be requested at info@dentaldirekt.de. We are continuously developing and enhancing our devices, and therefore reserve the right to make changes.

You will achieve the best color match by individually positioning the restoration in the blank (match). Depending on the height of the restoration, the positioning of the incisal, intermediate layer and dentin can be individually adjusted in the DD smart CAM 2.0 software to achieve the best possible color gradient.

This version replaces all previous versions.

Example of optimum nesting:

Indication: Crown, tooth 46 (lower first molar), shade A3

Blank: Zirkon BioStar Ultra ML, height 14 mm, color A3

CAM Software: DD smart CAM 2.0

Incisal	Blank position 1,0 mm		Blank position 2,5 mm		Blank position 3,3 mm	
	Transitional layers		Dentine		Dentine	
Dentine	More intense brightness in the body of the crown also. Shade is closer to A2.		Brightening in the cusps. Color gradient perfectly matched to the ratio of crown.		No brightening in the cusps. Chroma is darkened.	
	too light		match		too dark	

8.2 Milling, sintering and finishing

The blanks must only be processed with the milling systems intended for this purpose. The specifications of the machine manufacturer must be observed.

For the highest accuracy of fit, the specific magnification factor is printed on the side of the blank (disc) as the code to be used.

After the milling process, the frameworks must be checked for any visual defects. Damaged or contaminated frameworks must not be processed any further. The white frameworks can be colored with DD Basic Shade, DD Pro Shade C, DD Pro Shade Z and DD Art Elements effect colors, and constructions made of pre-colored zirconium oxide can be customized with DD Art Elements before sintering to full density (Manufacturer: Dental Direkt GmbH; observe separate instructions for use).

7. Handling and storage
Before using the material for the first time, check that the packaging and the blank itself are intact. Check whether the contents of the packaging correspond to the declaration on the label. Never use damaged material. Store the blanks only in the original packaging and in a cool, dry place. Avoid vibration, contamination and all contact with fluids.

8. Instructions for use in the laboratory
8.1 Processing / construction
Slident zirconium dioxide is a sensitive, high-performance ceramic and should be processed with special care, including in the partially sintered state!

The following construction parameters must always be considered when fabricating zirconium dioxide constructions:

Indication	Zirkon BioStar (colour), Zirkon BioStar Ultra (colour, ML)		Zirkon BioStar HT Smile (colour, ML)	
	Minimum wall thickness [mm]	Connector cross section [mm ²]	Minimum wall thickness [mm]	Connector cross section [mm ²]
Single crown	okkusal	0,5	0,5	–
	okkusal	0,5	0,5	–
Telescope	okkusal	0,7	0,7	–
	okkusal	0,5	0,5	–
Anterior Bridge 3 units	okkusal	0,5	0,8	> 10
	okkusal	0,5	0,8	> 12
Posterior Bridge 3 units	okkusal	0,7	0,8	–
	okkusal	0,7	0,8	> 12
Anterior Bridge 4 or more units	okkusal	0,6	0,8	> 9
	okkusal	0,8	0,8	> 12
Posterior Bridge 4 or more units	okkusal	0,7	0,8	–
	okkusal	0,8	0,8	> 12
Cantilever Bridge with 1 Pontic	okkusal	0,7	0,8	> 12
	okkusal	0,7	0,8	> 12

Weight per unit [g]

< 1	2	3	> 4
-----	---	---	-----

Heating and cooling rate [°C/min]

55	45	35	25
----	----	----	----

9. Instructions for use in the dental practice

For luting, we recommend conventional cementing with zinc oxide phosphate cements or glass ionomer cements. Luting composites can also be used. Ensure sufficient retention and a minimum stump height of 3 mm. For additional cleaning, the inner surface being bonded may be blasted with aluminum oxide (50 µm at 1-2 bar). Ensure that the surface is free of grease. Temporary luting is not recommended!

10. Material

	Zirkon BioStar (colour)	Zirkon BioStar Ultra (colour, ML)	Zirkon BioStar HT Smile (colour, ML)
ZrO ₂ + HfO ₂ + Y ₂ O ₃	≥ 99,0	≥ 99,0	≥ 99,0
Y ₂ O ₃	< 6	< 8	≤ 10
Al ₂ O ₃	< 0,5	< 0,15	≤ 0,01
Other oxides	< 1	< 1	< 1

Physical properties

	Zirkon BioStar (colour)	Zirkon BioStar Ultra (colour, ML)	Zirkon BioStar HT Smile (colour, ML)
CTE (25-500°C)	[10-6 K-1]	–10,8	–10,2
Chem. solubility	[µm/cm ²]	≤ 9,1	≤ 2,9
Fracture toughness (KIC)	[MPa·m]	> 9,5	4,0
Biegefestigkeit*	[MPa]	1150 ± 200	1200 ± 150

* measured according to DIN EN ISO 6872

11. Possible side effects and interactions
No known side effects or interactions.

12. Disposal
In compliance with local regulations. Non-contaminated and completely emptied packaging can be recycled.

Please note: Observe the information in the latest version of the safety data sheet.

13. Reporting incidents
Any serious incident that occurs in connection with the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

Please note: The safety and clinical performance summary report can be requested at info@dentaldirekt.de. We are continuously developing and enhancing our devices, and therefore reserve the right to make changes.

You will achieve the best color match by individually positioning the restoration in the blank (match). Depending on the height of the restoration, the positioning of the incisal, intermediate layer and dentin can be individually adjusted in the DD smart CAM 2.0 software to achieve the best possible color gradient.

This version replaces all previous versions.

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SILADENT Dr. Eberhard & Schöps GmbH