INSTRUCTIONS FOR USE

Dental LT Comfort Resin
V1 / V1.1

FLDLC001 / FLDLC011

EN	Instructions for Use	NO	Instruksjoner for bruk
DE	Gebrauchsanweisung	RO	Instrucțiuni de utilizare
FR	Mode d'emploi	TR	Kullanım talimatları
IT	Istruzioni per l'uso	LV	Lietošanas instrukcija
ES	Instrucciones de uso	LT	Naudojimo instrukcijos
PT	Instruções de uso	ET	Kasutamise juhised
NL	Gebruiksaanwijzing	SK	Návod na použitie
DA	Brugsanvisning	SL	Navodila za uporabo
SV	Bruksanvisning	PL	Instrukcja użytkowania
FI	Käyttöohjeet	BG	Инструкции за употреба
CS	Návod k použití	KO	사용 지침
HU	Használati útmutató	JA	使用説明書
EL	Οδηγίες χρήσης	ZH	使用说明

Prepared: 10/29/2024 PRNT-0080 Rev 05

INSTRUCTIONS FOR USE

Dental LT Comfort Resin
V1 / V1.1

FLDLC001 / FLDLC011

Keep away from sunlight

Use-by date

Consult Instructions for Use

European Conformity

Batch Code

Catalog Number

Manufacturer

Temperature Limit

Prescription only

Medical Device

Environmental hazard

Warning: Read Safety Information

Unique device identifier

Importer

Indicates the authorized
representative in Switzerland

Authorized representative
in the European Community

Formlabs Ohio Inc.
27800 Lemoyne Rd,
Millbury, OH 43447
USA
+1 617 855 0762

EU Importer
Formlabs GmbH
Mühlenstraße 15, Floor 1
10243 Berlin, Germany
+49 3091 734306

Curaden AG
Amlehnstrasse 22
6010 Kriens
Switzerland

Advena Ltd.
Tower Business Centre
2nd Flr. Tower Street,
Swatar, BKR 4013 Malta

Dental LT Comfort Resin V1.1 FLDLC011 can only be used with Form 4 generation printers.
 Dental LT Comfort Resin V1 FLDLC001 can only be used with Form 3 generation printers.
 Before printing, refer to the Manufacturing Guide and dental applications guide for this material at support.formlabs.com for printing and post-processing recommendations and requirements to ensure the correct and safe usage of this material. For detailed safety and environmental information, the Safety Data Sheet is available at support.formlabs.com.

I. Introduction And Indications For Use

INTRODUCTION AND INTENDED PURPOSE

Dental LT Comfort Resin is a light-curable polymer-based resin designed for the fabrication of biocompatible, long-term use, removable dental and orthodontic appliances, such as occlusal splints, nightguards, mouthguards, and/or repositioners by additive manufacturing.

INDICATIONS FOR USE

Dental LT Comfort Resin when utilized to print dental or orthodontic appliances such as occlusal splints, nightguards, mouthguards or repositioners is indicated to treat patients with conditions associated with occlusal wear such as Temporomandibular Joint Disorders (TMD) and/or Bruxism, respectively.

CONTRAINDICATIONS

Dental LT Comfort Resin should not be used in case of any known allergies or hypersensitivity history to photopolymer resin materials. In case of an allergic reaction, stop use immediately, and consult a physician.

INTENDED PATIENT POPULATION

All patients and populations in which a clinician has prescribed a treatment associated with the intended purpose.

INTENDED USERS

Dental and orthodontic professionals who fabricate dental and orthodontic appliances.

CLINICAL BENEFIT

The clinical benefit of Dental LT Comfort Resin is the ability to 3D print dental or orthodontic appliances such as occlusal splints, nightguards, mouthguards, and repositioners.

II. Safety

RESIN ASSOCIATED WARNINGS

Dental LT Comfort Resin (uncured) contains polymerizable monomers.

- Please read the Dental LT Comfort Resin SDS at support.formlabs.com for hazards associated with liquid resin.

RESIN USE ASSOCIATED PRECAUTIONS

- Protective glasses and nitrile gloves should be worn while handling resin. For information about handling Dental LT Comfort Resin, please read the SDS at support.formlabs.com.

INFORMATION FOR PATIENT SAFETY

- The patient should inform the clinician of any known allergy or hypersensitivity history to photopolymer resins.
- In case of an allergic reaction, stop use immediately, and consult a medical physician.
- Do not use appliances if damage or cracks are present.
- If signs of extreme wear, deterioration, and/or perforations appear, the appliances should be replaced to avoid material failure and/or potential bodily injury.
- Device lifetime will vary from patient to patient, however, devices should be replaced at least every six months. To maximize the device's lifetime, patients should be instructed to handle devices carefully and follow the cleaning and storage instructions.

III. Performance Characteristics

For detailed information on the mechanical properties of this material, technical data sheets can be found at [dental.formlabs.com](https://support.formlabs.com). Dental LT Comfort Resin has been evaluated in accordance with ISO 10993-1:2020, and ISO 7405:2018 and meets the requirements of all applicable biocompatibility testing.

IV. Specific Manufacturing Considerations

Refer to the Manufacturing Guide at support.formlabs.com for printing recommendations and requirements to ensure the correct and safe usage of this material.

- Hardware: Formlabs 3D Printer and accessories
- Software: Formlabs PreForm
- Part Orientation: Design guidelines for part orientation during printing can be found in the Manufacturing Guide.
- Part Thickness: Design guidelines for part thickness can be found in Manufacturing Guide
- Recommended Post-Processing Equipment:
 - Formlabs validated accessories
 - Formlabs validated washing units with approved solvents
 - Formlabs validated curing units

REQUIREMENTS

For biocompatibility compliance, Dental LT Comfort Resin was validated with a dedicated resin tank and mixer, build platform, wash unit, and post-processing equipment that were not mixed with any other resins.

A. PRINTING

Printing guidelines and validated equipment can be found in the Manufacturing Guide at support.formlabs.com.

B. WASHING

Washing guidelines and validated equipment can be found in the Manufacturing Guide.

C. POST-CURING

Post-Curing guidelines and validated equipment can be found in the Manufacturing Guide.

D. SUPPORT REMOVAL

Support marks can cause abrasion if not removed and polished. See the Manufacturing Guide for polishing recommendations.

E. CLEANING & DISINFECTION

1. Appliances may be cleaned using neutral soap and water or effervescent dental appliance cleaning tablets (used according to the manufacturer's directions).
2. Appliances may be disinfected by soaking in 70% IPA for five minutes per FDA guidelines.
3. Inspect appliances for cracks after cleaning/disinfection. Discard if damage or cracks are detected.

F. STORAGE

1. Printed appliances should be stored in closed, opaque or amber containers. Excess light exposure over time may affect the color and performance of printed appliances.
2. Store cartridges at 10°C - 25°C (50°F - 77°F). Do not exceed 25°C (77°F) when in storage.
3. Keep cartridges closed and away from ignition sources.

G. DISPOSAL

1. Cured resin is non-hazardous and may be disposed of as regular waste.
2. Follow facility protocols for waste that may be considered biohazardous.
3. Unused liquid resin & solvents should be disposed of in accordance with local regulations.
4. Cartridge & contaminated packaging should be disposed of in accordance with local regulations.
 - a. Contact a licensed professional waste disposal service to dispose of liquid resin.
 - b. Do not allow waste to enter storm or sewer drainage systems.
 - c. Avoid release into the environment.
 - d. Contaminated packaging: Dispose of cartridge or any unused product per local waste management guidelines.

H. ADVERSE EVENTS

To report any serious injury (temporary or permanent serious deterioration of health, a serious public health threat, or death), contact Formlabs Inc. (+1 617 855 0762) or your local Competent Authority.

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