
INSTRUCTIONS FOR USE

Permanent Crown Resin

Resin for 3D printing of permanent single crowns, inlays, onlays and veneers

|    |                      |    |                        |
|----|----------------------|----|------------------------|
| EN | Instructions for Use | LT | Naudojimo instrukcijos |
| DE | Gebrauchsanweisung   | SK | Návod na pouitie       |
| FR | Mode d'emploi        | SL | Navodila za uporabo    |
| IT | Istruzioni per l'uso | PL | Instrukcja użytkowania |
| ES | Instrucciones de uso | BG | Инструкции за употреба |
| PT | Instruções de uso    |    |                        |
| NL | Gebruiksaanwijzing   |    |                        |
| SV | Bruksanvisning       |    |                        |
| CS | Návod k použití      |    |                        |
| HU | Használati útmutató  |    |                        |
| EL | Οδηγίες χρήσης       |    |                        |

Prepared: 14/01/2025 PRNT-106964 REF 20811/02 Rev 06

INSTRUCTIONS FOR USE


# Permanent Crown Resin

Resin for 3D printing of permanent single crowns, inlays, onlays and veneers

|                                                                                   |                              |                                                                                   |                       |
|-----------------------------------------------------------------------------------|------------------------------|-----------------------------------------------------------------------------------|-----------------------|
|  | Keep away from sunlight      |  | Use-by date           |
|  | Consult Instructions for Use |  | Caution               |
|  | Batch Code                   |  | Catalog Number        |
|  | Manufacturer                 |  | Temperature Limit     |
|  | European Conformity          |  | Medical Device        |
|  | Date of Manufacture          |  | Professional use only |

 BEGO Bremer Goldschlägerei  
Wilh. Herbst GmbH & Co. KG  
Wilhelm-Herbst-Str. 1  
28359 Bremen, Germany  
Tel. +49.421.20.28-0  
Fax +49.421.20.28-100

**Manufactured for:** Formlabs Inc.  
35 Medford Street,  
Suite #1,  
Somerville, MA  
Tel. +1.617.702.8483

## 1. Intended Use

Permanent Crown is a light-curing, free-flowing plastic based on methacrylic acid esters for the production of permanent single crowns, inlays, onlays and veneers.

## 2. Indication

Permanent Crown is indicated for the following:

- Tooth defects due to trauma, tooth decay, attrition, abrasion or erosion / bio corrosion
- Weakening and / or incomplete application of the tooth structure (genetic dysplasia)
- Bite position changes, shape and position changes

## 3. Contraindications

Known allergy to one or more ingredients. In cases of doubt, the allergy should be clarified and ruled out based on a specific test prior to the application of this product. Permanent Crown should not be used for purposes other than the production of permanent single crowns, inlays, onlays and veneers. Any deviation from these instructions for use can have negative effects on the chemical and physical quality of plastics made from Permanent Crown.

## 4. Safety instructions

Permanent Crown is produced and tested according to the most stringent quality standards. It may only be used by a dental professional. In order to ensure optimum further processing, please read the information contained in the instructions for use carefully. The improper use of Permanent Crown and failure to follow information can have a detrimental effect on the quality of objects produced from Permanent Crown. Nitrile gloves, safety goggles and a coat must be worn as a means of protection when handling the resin and the plastic that has not been post-cured. Conventional medical gloves do not offer any lasting protection against the sensitizing effect of methacrylates. If the product comes into contact with the glove, take the glove off and discard it, wash your hands immediately with water and soap and put on a new glove. Consult a physician in the event of an allergic reaction. The safety and care instructions set down in the Permanent Crown instructions for use and safety data sheet shall apply to the handling of liquid resin and printed objects that have not been post-cured (objects in the "green condition"). A dust mask must be worn too due to potential dust formation while the printed objects are being processed.

It is prohibited to use plastic parts made of Permanent Crown as auxiliary equipment for food and drinks applications.

## 5. Side effects and precautions

### PRECAUTIONS / PROTECTION

It is essential that protective clothing be worn when handling Permanent Crown. Safety goggles and nitrile gloves must be used. However, we cannot completely rule out the possibility of personal reactions to individual components in isolated cases. In such cases, the respective user should discontinue use of Permanent Crown. If intolerances or allergic reactions occur when it comes into contact with the patient, discontinue use of the material.

#### Contains:

WARNING

Esterification products of 4,4'-isopropylidiphenol, ethoxylated and 2-methylprop-2-enoic acid. Silanized dental glass, methyl benzoylformate, diphenyl(2,4,6-trimethylbenzoyl) phosphine oxide. Total content of inorganic fillers (particle size 0.7 µm) is 30–50% by mass.

#### Hazard statements as per SDS

- Causes skin irritation.
- May cause an allergic skin reaction.
- Causes serious eye irritation.
- May cause respiratory irritation.
- May cause long harmful effects to aquatic life.

#### Precautionary statements as per SDS

- Avoid breathing mist / vapours / spray.
- Avoid release to the environment.
- Wear protective gloves / protective clothing / eye protection / face protection.
- Call a POISON CENTER / doctor if you feel unwell.
- If skin irritation or rash occurs: Get medical advice / attention.
- If eye irritation persists: Get medical advice / attention.
- Dispose of contents / container as per local and national regulations

#### NOTICE:

Any serious incident that has occurred in relation to the device due to a malfunction should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

## 6. General information on handling

Permanent Crown is supplied in colours according to the VITA<sup>1</sup> classical shades.

### STORAGE

Permanent Crown must be stored in the original cartridge at room temperature approx. 22 °C in a dark, dry place. It must be ensured that the temperature does not drop below +4 °C and does not exceed +28 °C! The minimum shelf life date printed on the product must be observed. Do not use resin cartridges past the expiration date on the cartridge label. Perfect processing cannot be guaranteed if materials which have exceeded their minimum shelf life date are used. The completely cured print objects must be stored at room temperature and protected from sources of bright light.

## 7. Processing

Permanent Crown is one of the system components in the Formlabs 3D print system and has been optimized for use in the SLA 3D printers – Form 3B and Form 3B+.

Create the object (STL-file) using a commercial CAD software, which is intended for dental applications. When designing, observe the requirements for minimum wall thicknesses for finished restorations. The design must match the anatomical dimensions of the tooth to be restored. The ratio of the crown height to the height of the bonding surface of the abutment must not exceed a value of 1.6.

---

<sup>1</sup> This symbol is a registered trademark of a company not affiliated with Formlabs Inc.

**NOTICE:****SINGLE CROWNS, INLAYS, ONLAYS AND VENEERS**

|                                          |        |
|------------------------------------------|--------|
| Minimum wall thicknesses anterior teeth  | 1.0 mm |
| Minimal wall thicknesses posterior teeth | 1.0 mm |
| Minimal wall thicknesses cervical        | 1.0 mm |
| Minimal wall thickness veneer shell      | 0.5 mm |

All accessories used must be designated for Permanent Crown only. For full compliance and biocompatibility of restorations from Permanent Crown, the process requires a dedicated resin tank, Stainless Steel Build Platform or Build Platform 2, Form Wash, and finishing kit which should not be mixed with any other resins.

The ideal working temperature of Permanent Crown is in the temperature range between 18 °C and 28 °C . **Before using a new cartridge of Permanent Crown resin, the material has to be shaken well!** For further processing – selecting the resin, setting up the print job – as part of the printing process, follow the Formlabs Form 3B or Form 3B+ printer instructions for use. Permanent Crown in the resin tank does not need to be mixed prior to printing, up to a month of non-use, as solids (fillers) in the material do not settle.

**PRINTING:**

- Insert resin cartridge, Stainless Steel Build Platform or Build Platform 2, and the compatible resin tank into the Formlabs 3D printer.
- Prepare a print job using PreForm software. Import a dental restoration STL file. Orient files horizontally with the occlusal plane facing the build platform and generate supports. For recommendations on print orientation and support placement, see the detailed application guide at [support.formlabs.com](https://support.formlabs.com).
- Send print job to printer. Begin print by selecting print job from print menu. Follow any prompts or dialogs shown on the printer screen. Printer will automatically complete the print.

**8. Subsequent processing**

On completion of printing, the print objects are released from the build platform using the spatula supplied. The print objects should be cleaned using the Form Wash as follows.

**WASHING:**

- Remove the build platform with the printed parts from the printer and place it into a Form Wash with clean IPA ( $\geq 99\%$ ).
- Set the Form Wash for 3 minutes. Make sure the parts are fully submerged in IPA when washing. Exceeding wash duration may affect dimensional accuracy and performance of printed parts over time.
- Once parts are washed, inspect parts for any remaining uncured resin on print surfaces. Use an IPA-filled spray bottle to remove any remaining resin on the printed parts and in between supports and rafts. An IPA-soaked brush may be used to help with uncured resin removal.

**DRYING:** Use compressed air to dry the parts. The dried parts will have a white, powdery coating on the print surfaces. If any wet, uncured resin is still present after drying, reuse the IPA-soaked brush and IPA-filled squeeze bottle to remove uncured resin and dry again.

**REMOVAL:** After drying, remove printed parts from the build platform by wedging the part removal tool under the print raft and rotating the tool. For detailed techniques, see the support material on [dental.formlabs.com](https://dental.formlabs.com).

**POST CURE:** The final properties of the print object depend on the post-curing process to maintain dimensional accuracy and biocompatibility:

## 8.1. POST-CURING WITH FORM CURE

8.1.1. With the raft and supports still intact, place printed parts in the Form Cure with the raft side down. Form Cure parts at 60°C for 20 minutes.

- a. Remove parts from Form Cure and use a handpiece with a cutting disc or side cutters to separate supports and raft from printed parts.
- b. Carefully sandblast the printed part surfaces to remove the white coating using a sandblaster with glass bead blasting material 50µm, ( e.g. Perlablast® micro<sup>2</sup>, BEGO, REF 46092/54302) at a maximum blasting pressure of 1.5 bar. If a sandblaster is unavailable, refer to Formlabs application guide for an alternative manual cleaning method.
- c. Check printed part fit on the model and finish and contour as needed using a dental handpiece and carbide rotary burs or diamond grinding stones.

8.1.2. Place finished parts in the Form Cure and post cure again at 60°C for an additional 20 minutes, with the occlusal plane facing upwards. For best results and full compliance, printed parts must be post cured twice as recommended. Any deviations from this protocol, may lead to unsatisfactory fit and affect biocompatibility.

## 8.2. POST-CURING WITH FAST CURE

8.2.1. With the raft and supports still intact, place printed parts in the Fast Cure with the raft side down. Cure parts for 2.5 minutes at Light Intensity 1.

- a. Follow steps 8.1.1.a-8.1.1.c to remove supports and sandblast parts.

8.2.2. Place finished parts in the Fast Cure and post cure again at Light Intensity 1 for additional 2.5 minutes, with the occlusal plane facing upwards. For best results and full compliance, printed parts must be post cured twice as recommended. Any deviations from this protocol may lead to unsatisfactory fit and affect biocompatibility.

## 8.3. POST-CURING WITH FORM CURE (2ND GENERATION)

8.3.1. Follow steps 8.1.1.a-8.1.1.c to remove supports and sandblast parts.

8.3.2. Place finished parts in the Form Cure (2nd Generation) with the occlusal plane facing upwards. Post cure for 10 minutes at 60°C. Only one curing cycle is required with Form Cure (2nd Generation).

### NOTICE:

The times given only apply to regularly maintained equipment that guarantees a corresponding light intensity. Inspect the dental restorations for cracks. Discard if any cracks or damage are detected. If the post-curing process is interrupted by a malfunction, the printed object should not be used until it has been cured for one full cycle. Refer to the operating instructions for the light-curing device to find out how the error can be resolved, and then repeat the post-curing process with the printed objects.

### SUPPLEMENTING OF PRINTED OBJECTS (OPTIONAL)

Defects (e. g. missing contact points) can be supplemented with the resin or with commercially available composite veneering materials.

#### Supplementing with resin

- Blast the areas to be supplemented with aluminum oxide – pressure 1.5 bar and particle size 110 µm.
- Put some liquid Permanent Crown on the object.
- Put it under light for a short time (e.g. couple of minutes in the Form Cure).
- If more material has to be applied, put additional liquid Permanent Crown on the last layer and then put it again under light for a short time (e.g. couple of minutes in the Form Cure).
- Final curing of the objects (e.g. 2 x 20 minutes at 60°C in the Form Cure).

#### Supplementing with veneering material

The printed objects can be supplemented with commercially available composite veneering materials (e. g. VITAVM LC<sup>2</sup>, Vita Zahnfabrik, und VITAVM LC flow<sup>2</sup>, Vita Zahnfabrik). The manufacturer's instructions for use must be observed.

**CAUTION:**

The dental object may only be repaired or supplemented outside the patient's mouth and by a dental professional.

**POLISHING:** Polish the surface of the objects with pumice stone and polishing compound

- a. Avoid overheating the parts during polishing.
- b. Optimal surface quality is achieved by polishing after post-curing.
- c. If rough marks are left, sand and polish to smooth surfaces to improve patient comfort.

**TIP:**

Optionally, the surface of the objects can be coated with light-curing glaze (e. g. Vita ENAMIC GLAZE<sup>3</sup>, Vita Zahnfabrik or GC OPTIGLAZE<sup>3</sup>, GC). Pay attention to the manufacturer's instructions for use.

**INDIVIDUALIZATION (OPTIONAL)**

Instead of polishing the surface, the post-cured objects can also be individualized with composite stains or veneering materials.

**Individualization with composite stain**

Commercially available composite stain materials (e. g. VITA ENAMIC STAIN<sup>3</sup>, Vita Zahnfabrik and GC OPTIGLAZE Color<sup>3</sup>, GC) can be used for individualization of the post-cured Permanent Crown restoration. The manufacturer's instructions for use must be observed.

**Individualization with veneering materials**

Commercially available composite veneering materials (e. g. VITAVM LC<sup>3</sup>, Vita Zahnfabrik and VITAVM LC flow<sup>3</sup>, Vita Zahnfabrik) can be used for individualization of the restoration. The manufacturer's instructions for use must be observed.

**VENEER SHELLS**

Commercially available light-curing veneering composites can be used to attach the post-cured veneer shells to the metal framework (e. g. VITAVM LC, VITAVM LC flow, Vita Zahnfabrik).

**Basic work steps**

- Condition the surface of the metal framework.
- Cover the metal framework with opaquer.
- Sandblast the inside of the veneer shell and clean the veneer shell from dust.
- Inner surface conditioning of the veneer shell. Observe the information of the veneering material manufacturer in section "Individualization of composite frameworks" or "Individualization of artificial teeth".
- Applying the veneering material to the inside of the veneer shell or the metal framework.
- Positioning of the veneer shell on the metal framework and remove the excess with a brush or an instrument.
- Carry out the polymerization according to the instructions of the veneering material manufacturer.
- Finally, the surface of the veneer shell should be polished or customized with composite stains (see section "Polishing" or "Individualization").

**9. Cleaning post-cured objects (optional)**

Fully post-cured objects made from Permanent Crown can be steam cleaned or immersed in a bath of cleaning solution.

**10. Luting****CEMENTATION ON ABUTMENTS**

The definitive restoration can be fixed onto titanium abutments (e. g. BEGO Semados® solid titanium abutments) with Panavia v5 (Kuraray Noritake). The instructions for use of the luting agent manufacturer must be observed.

## ATTACHMENT TO TOOTH STUMP

The finished permanent restorations can be attached using self-adhesive cements (e. g. RelyX Unicem<sup>3</sup>, 3M Espe) or composite cement with a primer (e. g. Variolink Esthetic DC<sup>3</sup> and Monobond Plus<sup>\*</sup>, Ivoclar Vivadent).

### NOTICE:

- Observe the instructions for use of the luting agent. It is not required to etch the restoration before attaching.
- Additional exposure to curing lights after attachment will not affect the properties of the finished object.
- Crowns made of Permanent Crown are not suitable for use on titanium adhesive bases.

## 11. Disposal

The cured, separated material (base plate, support structure) can no longer be used. Cured material can be disposed of as domestic waste. Unused resin, contaminated packaging or IPA used for cleaning with resin residues must be disposed of via the local waste disposal authority or a hazardous waste collection point stating the safety data sheet.

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