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FotoDent® IBT

Intended use

Material for the manufacturing of dental indirect bonding trays.

Indication for use: Light-curing material for the manufacturing of dental indirect bonding trays using a 3D printing process for placing brackets.

Indication: For the treatment of misaligned teeth.

Contraindication: Intended use of FotoDent® IBT is for manufacturing dental bonding trays only; an allergic reaction might occur in case of any allergy to (meth)acrylates.

Patient target group: Patients with defective teeth positions.

Intended users: The dental technicians or specialist staff in the dental laboratory is responsible for the production of indirect bonding trays in 3D printing. The dentist or orthodontist uses the finished indirect bonding tray on the patient.

Requirements

Printer: Rapid Shape DII series, D20+, D30+ / Asiga MAX UV, PRO 4K Dreve: Phrozen Sonic XL 4K (qualified by Dreve) / Phrozen Sonic XL 4K 2022 (qualified by Dreve)

Software: Autodesk Netfabb®, Composer, Dreve ElementS

Post curing unit: PCU LED N₂, FotoDent® flash, Otofash G171, RS cure (see overview)

Basic Material Properties

Please refer to Technical Data Sheet (TDS). Available on request.

Specific Manufacturing Considerations

1. Digital denture indirect bonding tray file

1.1 File format: STL file

1.2 Digital design: dental indirect bonding tray

2. Asiga and software

2.1 Hardware: See directions for use in the printer manual.

2.2 Software: See directions for use in the printer manual.

2.3 Printing Parameters:

- a) Slice thickness: 100µm
- b) Optimal Orientation: horizontal, directly on the platform
- c) Support point size: no supports are needed

2.4 Environmental Conditions

- a) MAX UV print room temperature: 35 °C ± 3 °C
- b) Humidity: 30–70 %

3. Rapid Shape and software

3.1 Hardware: See directions for use in the printer manual.

3.2 Software: See directions for use in the printer manual.

3.3 Printing Parameters

- a) Slice thickness: 100 µm
- b) Optimal Orientation: horizontal, directly on the platform
- c) Support point size: no supports are needed

3.4 Environmental Conditions

- a) Print temperature: See directions for use in the printer manual.
- b) Humidity: 30–70 %

4. Dreve Phrozen Printer and software

4.1 Hardware: See directions for use in the printer manual. Firmware version: 1.62 or higher

Lot number / Expiry date

The lot number and the expiry date are indicated on each FotoDent® IBT packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

4.3 Printing Parameters

- a) Slice thickness: 100 µm
- b) Optimal Orientation: horizontal, directly on the platform
- c) Support point size: no supports are needed

4.4 Environmental Conditions

- a) Print temperature: See directions for use in the printer manual.
- b) Humidity: 30–70 %

5. Notification

The specifications and biocompatibility have been validated using printer software (Autodesk® Netfabb, Composer, Dreve ElementS) and aforementioned printer and process parameters indicated.

Caution

- Due to the flexibility of the FotoDent® IBT material it is recommended to build the object directly on the platform without support structures in an horizontal orientation. Therefore, the occlusal area of the indirect bonding tray has to be designed as a flat surface without chewing surface structure. Details can be found in the instruction for use of your design program.

- Due to different bracket shapes, please be aware of sufficient fitting during designing IBT. If not the brackets may fall out of lower jaw IBT during placing IBT.
- To avoid damage of material quality do not expose the liquid material to light under any circumstances. Variations from the described manufacturing process may impair biocompatibility, lead to modified mechanical characteristics and/ or color variations of the FotoDent® IBT material.

- Impurity due to the construction or a break of the material and impurity due to operation mistakes cannot be excluded. Thanks to the low viscosity however it is possible to filtrate FotoDent® IBT. It is recommended to regularly take out the container of the production unit, to homogenize and filtrate the content. Upcoming entrained bubbles get eliminated through a holding time of approx. 1–2 hours with FotoDent® IBT.

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