LIQCREATE

GB

LIQCREATE TRAY

Instructions for Use

General information
Liqcreate Tray is a methacrylate photoinitiator, pigments and additives based light-curing resin for use in the dental laboratory. For DLP/LCD 3D printers operating at a wavelength of 385-405nm.

Intended use
3D printing resin for the manufacture of individual impression trays.

Special precautions
- Non-compliance with instructions and improper use may lead to impairment of quality.
- Wear protective gloves, protective clothing and safety goggles while working.
- Avoid contact with eyes! If liquid material comes into contact with eyes inadvertently, rinse immediately with plenty of water and consult a doctor if necessary.
- Avoid direct skin contact with liquid material and components prior to post-curing.
- Keep out of reach of children.
- Observe safety data sheet.
- Report any serious incident that has occurred in relation to the device to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Side effect
- Irritation of eyes and skin.
- The polymerised material should not be used if there are allergies to any of the ingredients or in case of contact allergies.
- In rare cases, allergic reactions to components of Liqcreate Tray may occur.
- Inhalation of monomer vapour and skin contact with non-polymerised materials should be avoided.

Notes on use
- Shake the bottle for at least one minute before use.
- When filling the vat with resin, observe the instructions of the device manufacturer.
- Note the temperature that the device manufacturer specifies for production. The resin can also be processed at room temperature.
- Darken the room slightly before removing the print platform with the printed object. Then, open the hood and remove the printing platform with the printed object.
- Remove the printed object from the platform.
- Only remove the support structures after cleaning and post-curing.

Reusing used resin
Sieve the remaining resin in the vat first (removal of polymerised particles) with a sieve gauze (about 200 µm). After sieving, the used resin can be mixed with fresh resin 1:1 and reused in the printer.

Post-processing of the printed object
- Clean the printed object with ethanol and a soft brush for a maximum of 3 minutes.
- Optional: Clean printed object in a separate container with fresh ethanol in an ultrasonic bath as needed (5-10 mins).
- Optional: Clean openings, holes and gaps with compressed air.

Post-exposure of the printed object
It is recommended to post-cure the printed object on each side for 20 minutes in a UVA light-curing device with 4 lamps of 9 watts each. Please note the information provided by the device manufacturer.

Restrictions on use
- Intended exclusively for the creation of dental trays in the dental sector.
- Must not be brought into direct contact with the patient.

Classification
Class I medical product in accordance with Regulation (EU) 2017/745.

Composition
Mixture of acrylate and methacrylate resins, photoinitiators, pigments and additives.

Storage and durability
- Store at room temperature/0-22 °C in a dark, dry place.
- Only handle Liqcreate Tray as briefly as possible to reduce exposure to daylight/room light. Excessive exposure to light leads to an unwanted polymerization of the material.
- Close the container immediately after each use and keep it tightly closed.
- Do not use after the expiration date.

Packaging
REF LTR01000 Liqcreate Tray 1000g

Liqcreate
Texasdreef 7
3565 CL Utrecht
Netherlands
Tel: +(31) 085 060 58 49
www.liqcreate.com
info@liqcreate.com

dent a pharm
Produktionsgesellschaft mbH
Schusterring 35
25355 Barmstedt / Germany
Tel: +(49) 4123/9225-0
Fax: +(49) 4123/9228-49
www.dent-a-pharm.de
info@dent-a-pharm.de

NL

LIQCREATE TRAY

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Algemene informatie
Liqcreate Tray is a methacrylate photoinitiator, pigments and additives based light-curing resin for use in the dental laboratory. For DLP/LCD 3D printers operating at a wavelength of 385-405nm.

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